



**Product Specification**  
**Sarpette® M, single-channel, variable, 100 µl - 1,000 µl, blue, 1 piece(s)/case**



**Product description**

Order number 90.3100.000  
Product description Sarpette® M, single-channel, variable, 100 µl - 1,000 µl, blue, 1 piece(s)/case

**Product characteristics**

Volume of work 100 µl - 1.000 µl  
Coding colour blue  
Autoclavable (max. °C) 121 °C

**Material & colours**

Colour of product grey

**Purity & certification**

Product category no medical device | IVD  
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Sarpette® M, single-channel, variable, 100 µl -**  
**1,000 µl, blue, 1 piece(s)/case**

**Page 2**

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Depth of box	315 mm
Width of box	90 mm
Height of box	30 mm
Weight of product	0.088 kg
EAN of inner box	4038917937061
EAN case	4038917937078

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**