



Product description

Order number	90.3100.200
Product description	Sarpette® M, single-channel, variable, 20 µl - 200 µl, yellow, 1 piece(s)/case

Product characteristics

Volume of work	20 µl - 200 µl
Coding colour	yellow
Autoclavable (max. °C)	121 °C

Material & colours

Colour of product	grey
-------------------	------

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Depth of box	315 mm
Width of box	90 mm
Height of box	30 mm
Weight of product	0.088 kg
EAN of inner box	4038917937009
EAN case	4038917937016

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

