



Product description

| | |
|---------------------|--|
| Order number | 90.3112.010 |
| Product description | Sarpette® M, 12-channel, variable, 0.5 µl - 10 µl, grey, 1 piece(s)/case |

Product characteristics

| | |
|------------------------|----------------|
| Volume of work | 0,5 µl - 10 µl |
| Coding colour | grey |
| Autoclavable (max. °C) | 121 °C |

Material & colours

| | |
|-------------------|------|
| Colour of product | grey |
|-------------------|------|

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 1 |
| Piece(s) / pallet | 1 |
| Depth of box | 288 mm |
| Width of box | 150 mm |
| Height of box | 28 mm |
| Weight of product | 0.185 kg |
| EAN of inner box | 4038917937610 |
| EAN case | 4038917937627 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

