




**Product Specification**  
**Serological pipette, with tip, plugged, 1 ml, sterile,  
25 piece(s)/bag**

**Page 1**



**Product description**

Order number 86.1251.025  
Product description Serological pipette, with tip, plugged, 1 ml, coding colour:  
yellow, material: PS, sterile, 25 piece(s)/bag

**Product characteristics**

Label/ Print with printed graduations  
Colour of print/label yellow/black  
Graduation yes  
Design with tip, plugged  
Volume of work 0,1 ml - 1,3 ml  
Coding colour yellow

**Size**

Sample volume 1 ml  
Diameter 4 mm  
Length of product 270 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Serological pipette, with tip, plugged, 1 ml, sterile,**  
**25 piece(s)/bag**

**Page 2**

**Material & colours**

Product material Polystyrene (PS)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Purity standard sterile  
Sterilisation Gamma irradiation  
Batched yes

**Packaging**

Minimum order qty. 1000  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 25  
Piece (s) / outer case 1000  
Piece(s) / pallet 108000  
Depth of case 413 mm  
Width of case 250 mm  
Height of case 150 mm  
Case volume 0.0155 cbm  
Weight of product 0.0033 kg  
Weight of case 3.94 kg  
EAN of inner box 4038917127165  
EAN case 4038917021074

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**