




Product Specification
Serological pipette, with tip, plugged, 2 ml, sterile,
25 piece(s)/bag

Page 1



Product description

Order number	86.1252.025
Product description	Serological pipette, with tip, plugged, 2 ml, coding colour: green, material: PS, sterile, 25 piece(s)/bag

Product characteristics

Label/ Print	with printed graduations
Colour of print/label	green/black
Graduation	yes
Design	with tip, plugged
Volume of work	0,2 ml - 2,5 ml
Coding colour	green

Size

Sample volume	2 ml
Diameter	6 mm
Length of product	270 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Serological pipette, with tip, plugged, 2 ml, sterile,
25 piece(s)/bag

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard sterile
Sterilisation Gamma irradiation
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest bag
subpackaging
Piece(s) / inner box 25
Piece (s) / outer case 1000
Piece(s) / pallet 81000
Depth of case 413 mm
Width of case 250 mm
Height of case 195 mm
Case volume 0.0201 cbm
Weight of product 0.0036 kg
Weight of case 4.21 kg
EAN of inner box 4038917127219
EAN case 4038917021098

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT