




Product Specification
**Serological pipette, without tip, plugged, 5 ml,
sterile, 10 piece(s)/bag**

Page 1



Product description

Order number	86.1687.010
Product description	Serological pipette, without tip, plugged, 5 ml, material: PS, sterile, 10 piece(s)/bag

Product characteristics

Label/ Print	with printed graduations
Colour of print/label	black
Graduation	yes
Design	without tip, plugged
Volume of work	0,6 ml - 5 ml

Size

Sample volume	5 ml
Diameter	8 mm
Length of product	340 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Serological pipette, without tip, plugged, 5 ml,
sterile, 10 piece(s)/bag**

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard sterile
Sterilisation Gamma irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest bag
subpackaging
Piece(s) / inner box 10
Piece (s) / outer case 500
Piece(s) / pallet 40500
Depth of case 413 mm
Width of case 250 mm
Height of case 195 mm
Case volume 0.0201 cbm
Weight of product 0.0064 kg
Weight of case 3.88 kg
EAN of inner box 4038917133074
EAN case 4038917024686

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT