



### **Product description**

Order number	86.1170
Product description	Transfer pipette, nominal volume: 1 ml, (LxW): 115 x 10 mm, material: LD-PE, transparent, droplet size: 40 µl - 50 µl, injected graduation up to 0.3 ml, 500 piece(s)/case

### **Product characteristics**

Graduation	yes
------------	-----

#### **Size**

Sample volume	1 ml
Width of product	10 mm
Length of product	115 mm
Drop size	40 µl - 50 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 1 ml, (LxW): 115 x 10 mm, LD-PE,  
transparent**

**Page 2**

**Material & colours**

Product material Low Density Polyethylene (LD-PE)  
Colour of product transparent

**Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

**Packaging**

Minimum order qty. 2000  
Type of smallest case  
subpackaging  
Piece(s) / inner box 500  
Piece (s) / outer case 2000  
Piece(s) / pallet 168000  
Depth of box 261 mm  
Width of box 126 mm  
Height of box 121 mm  
Depth of case 533 mm  
Width of case 263 mm  
Height of case 140 mm  
Case volume 0.0196 cbm  
Weight of product 0.0007 kg  
Weight of case 1.89 kg  
EAN of inner box 4038917085564  
EAN case 4038917004206

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**