



Product description

Order number	86.1170.300
Product description	Transfer pipette, nominal volume: 1 ml, (LxW): 115 x 10 mm, material: LD-PE, transparent, droplet size: 40 µl - 50 µl, 500 piece(s)/case

Product characteristics

Graduation	yes
------------	-----

Size

Sample volume	1 ml
Width of product	10 mm
Length of product	115 mm
Drop size	40 µl - 50 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Transfer pipette, 1 ml, (LxW): 115 x 10 mm, LD-PE,
transparent**

Page 2

Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 500
Type of smallest case
subpackaging
Piece(s) / inner box 500
Piece (s) / outer case 4000
Piece(s) / pallet 156000
Case volume 0.0401 cbm
Weight of product 0.001 kg
Weight of case 4.55 kg
EAN of inner box 4038917378352
EAN case 4038917378345

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT