




**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,**  
**LD-PE, transparent**

**Page 1**



**Product description**

Order number	86.1171
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 155 x 15 mm, material: LD-PE, transparent, droplet size: 35 µl - 55 µl, injected graduation up to 3 ml, 500 piece(s)/case

**Product characteristics**

Graduation	yes
------------	-----

**Size**

Sample volume	3.5 ml
Width of product	15 mm
Length of product	155 mm
Drop size	35 µl - 55 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,**  
**LD-PE, transparent**

**Page 2**

**Material & colours**

Product material Low Density Polyethylene (LD-PE)  
Colour of product transparent

**Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

**Packaging**

Minimum order qty. 2000  
Type of smallest case  
subpackaging  
Piece(s) / inner box 500  
Piece (s) / outer case 2000  
Piece(s) / pallet 56000  
Depth of box 261 mm  
Width of box 161 mm  
Height of box 231 mm  
Depth of case 533 mm  
Width of case 333 mm  
Height of case 250 mm  
Case volume 0.0444 cbm  
Weight of product 0.0015 kg  
Weight of case 3.92 kg  
EAN of inner box 4038917085571  
EAN case 4038917004213

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-09-26

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**