




Product Specification
Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm,
LD-PE, transparent

Page 1



Product description

Order number	86.1172
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 156 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl, injected graduation up to 1 ml, 500 piece(s)/case

Product characteristics

Graduation	yes
------------	-----

Size

Sample volume	3.5 ml
Width of product	12.5 mm
Length of product	156 mm
Drop size	30 µl - 45 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm,
LD-PE, transparent

Page 2

Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 2000
Type of smallest case
subpackaging
Piece(s) / inner box 500
Piece (s) / outer case 2000
Piece(s) / pallet 56000
Depth of box 261 mm
Width of box 161 mm
Height of box 231 mm
Depth of case 533 mm
Width of case 333 mm
Height of case 250 mm
Case volume 0.0444 cbm
Weight of product 0.0013 kg
Weight of case 3.47 kg
EAN of inner box 4038917085588
EAN case 4038917004220

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT