




Product Specification
Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm,
LD-PE, transparent

Page 1



Product description

Order number	86.1172.010
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 156 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl, injected graduation up to 1 ml, sterile, 10 piece(s)/blister

Product characteristics

Graduation	yes
------------	-----

Size

Sample volume	3.5 ml
Width of product	12.5 mm
Length of product	156 mm
Drop size	30 µl - 45 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm,
LD-PE, transparent

Page 2

Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Gas sterilisation
Batched yes

Packaging

Minimum order qty. 1200
Type of smallest blister
subpackaging
Piece(s) / inner box 300
Piece (s) / outer case 1200
Piece(s) / pallet 33600
Depth of box 261 mm
Width of box 161 mm
Height of box 231 mm
Depth of case 533 mm
Width of case 333 mm
Height of case 250 mm
Case volume 0.0444 cbm
Weight of product 0.0013 kg
Weight of case 2.79 kg
EAN of inner box 4038917353151
EAN case 4038917353120

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT