



### **Product description**

Order number	86.1172.300
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 156 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl, 500 piece(s)/case

### **Product characteristics**

Graduation	yes
------------	-----

### **Size**

Sample volume	3.5 ml
Width of product	12.5 mm
Length of product	156 mm
Drop size	30 µl - 45 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm,**  
**LD-PE, transparent**

**Page 2**

**Material & colours**

Product material Low Density Polyethylene (LD-PE)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest case  
subpackaging  
Piece(s) / inner box 500  
Piece (s) / outer case 5000  
Piece(s) / pallet 56000  
Case volume 0.0994 cbm  
Weight of product 0.0018 kg  
Weight of case 10 kg  
EAN of inner box 4038917378413  
EAN case 4038917378406

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**