




**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm,**  
**LD-PE, transparent**

**Page 1**



### Product description

Order number	86.1172.301
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 156 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl, sterile, 1 piece(s)/blister

### Product characteristics

Graduation	yes
------------	-----

### Size

Sample volume	3.5 ml
Width of product	12.5 mm
Length of product	156 mm
Drop size	30 µl - 45 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm,**  
**LD-PE, transparent**

**Page 2**

**Material & colours**

Product material Low Density Polyethylene (LD-PE)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

**Packaging**

Minimum order qty. 4000  
Type of smallest blister  
subpackaging  
Piece(s) / inner box 500  
Piece (s) / outer case 4000  
Piece(s) / pallet 30240  
Case volume 0.1359 cbm  
Weight of product 0.0032 kg  
Weight of case 13.09 kg  
EAN of inner box 4038917378444  
EAN case 4038917378437

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**