



### **Product description**

Order number	86.1174
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 152 x 15 mm, material: LD-PE, transparent, droplet size: 40 µl - 55 µl, 400 piece(s)/case

### **Product characteristics**

Graduation	no
------------	----

### **Size**

Sample volume	6 ml
Width of product	15 mm
Length of product	152 mm
Drop size	40 µl - 55 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Low Density Polyethylene (LD-PE)  
Colour of product transparent

### **Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

### **Packaging**

Minimum order qty. 1600  
Type of smallest case  
subpackaging  
Piece(s) / inner box 400  
Piece (s) / outer case 1600  
Piece(s) / pallet 44800  
Depth of box 261 mm  
Width of box 161 mm  
Height of box 231 mm  
Depth of case 533 mm  
Width of case 333 mm  
Height of case 250 mm  
Case volume 0.0444 cbm  
Weight of product 0.0017 kg  
Weight of case 3.58 kg  
EAN of inner box 4038917085601  
EAN case 4038917004244

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

