



### **Product description**

Order number	86.1174.020
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 152 x 15 mm, material: LD-PE, transparent, droplet size: 40 µl - 55 µl, sterile, 20 piece(s)/bag

### **Product characteristics**

Graduation	no
------------	----

### **Size**

Sample volume	6 ml
Width of product	15 mm
Length of product	152 mm
Drop size	40 µl - 55 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 6 ml, (LxW): 152 x 15 mm, LD-PE,  
transparent**

**Page 2**

**Material & colours**

Product material Low Density Polyethylene (LD-PE)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

**Packaging**

Minimum order qty. 3200  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 400  
Piece (s) / outer case 3200  
Piece(s) / pallet 48000  
Case volume 0.1395 cbm  
Weight of product 0.0023 kg  
Weight of case 8.95 kg  
EAN of inner box 4038917318402  
EAN case 4038917318396

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**