



Product description

Order number	86.1175
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 146 x 15 mm, material: LD-PE, transparent, droplet size: 20 μ l - 30 μ l, 400 piece(s)/case

Product characteristics

Graduation	no
------------	----

Size

Sample volume	6 ml
Width of product	15 mm
Length of product	146 mm
Drop size	20 μ l - 30 μ l

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Transfer pipette, 6 ml, (LxW): 146 x 15 mm, LD-PE,
transparent**

Page 2

Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 1600
Type of smallest case
subpackaging
Piece(s) / inner box 400
Piece (s) / outer case 1600
Piece(s) / pallet 44800
Depth of box 261 mm
Width of box 161 mm
Height of box 231 mm
Depth of case 533 mm
Width of case 333 mm
Height of case 250 mm
Case volume 0.0444 cbm
Weight of product 0.0015 kg
Weight of case 3.33 kg
EAN of inner box 4038917085618
EAN case 4038917004251

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT