



### **Product description**

Order number	86.1175.001
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 146 x 15 mm, material: LD-PE, transparent, droplet size: 20 µl - 30 µl, sterile, 1 piece(s)/blister

### **Product characteristics**

Graduation	no
------------	----

### **Size**

Sample volume	6 ml
Width of product	15 mm
Length of product	146 mm
Drop size	20 µl - 30 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Low Density Polyethylene (LD-PE)  
Colour of product transparent

### **Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Gas sterilisation  
Batched yes

### **Packaging**

Minimum order qty. 840  
Type of smallest blister  
subpackaging  
Piece(s) / inner box 42  
Piece (s) / outer case 840  
Piece(s) / pallet 23520  
Depth of case 533 mm  
Width of case 333 mm  
Height of case 250 mm  
Case volume 0.0444 cbm  
Weight of product 0.0016 kg  
Weight of case 2.56 kg  
EAN of inner box 4038917126311  
EAN case 4038917020589

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

