




Product Specification
Transfer pipette, 2 ml, (LxW): 154 x 11 mm, LD-PE,
transparent

Page 1



Product description

Order number	86.1176.300
Product description	Transfer pipette, nominal volume: 2 ml, (LxW): 154 x 11 mm, material: LD-PE, transparent, droplet size: 35 µl - 45 µl, 500 piece(s)/case

Product characteristics

Graduation	no
------------	----

Size

Sample volume	2 ml
Width of product	11 mm
Length of product	154 mm
Drop size	35 µl - 45 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Transfer pipette, 2 ml, (LxW): 154 x 11 mm, LD-PE,
transparent**

Page 2

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	case
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	84000
Case volume	0.0994 cbm
Weight of product	0.0015 kg
Weight of case	8.64 kg
EAN of inner box	4038917378567
EAN case	4038917378550

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT