



### **Product description**

Order number	86.1180
Product description	Transfer pipette, nominal volume: 1 ml, (LxW): 104 x 10 mm, material: LD-PE, transparent, droplet size: 17 $\mu$ l - 22 $\mu$ l, 750 piece(s)/case

### **Product characteristics**

Graduation	no
------------	----

### **Size**

Sample volume	1 ml
Width of product	10 mm
Length of product	104 mm
Drop size	17 $\mu$ l - 22 $\mu$ l

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 1 ml, (LxW): 104 x 10 mm, LD-PE,  
transparent**

**Page 2**

**Material & colours**

Product material Low Density Polyethylene (LD-PE)  
Colour of product transparent

**Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

**Packaging**

Minimum order qty. 3000  
Type of smallest case  
subpackaging  
Piece(s) / inner box 750  
Piece (s) / outer case 3000  
Piece(s) / pallet 252000  
Depth of box 261 mm  
Width of box 126 mm  
Height of box 121 mm  
Depth of case 533 mm  
Width of case 263 mm  
Height of case 140 mm  
Case volume 0.0196 cbm  
Weight of product 0.0005 kg  
Weight of case 2.07 kg  
EAN of inner box 4038917085649  
EAN case 4038917004282

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**