



### **Product description**

Order number	93.1646.002
Product description	DishRack, height: 370 mm, material: PC, yellow, for 88 petri dishes with 92 mm Ø, 1 piece(s)/case

### **Product characteristics**

Storage capacity	88
------------------	----

### **Size**

Height of product	370 mm
-------------------	--------

### **Material & colours**

Product material	Polycarbonate (PC)
Colour of product	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category no medical device | IVD  
Batched no

**Packaging**

Minimum order qty. 1  
Type of smallest subpackaging case  
Piece(s) / inner box 1  
Piece (s) / outer case 1  
Piece(s) / pallet 108  
Depth of case 383 mm  
Width of case 193 mm  
Height of case 190 mm  
Case volume 0.014 cbm  
Weight of product 1.086 kg  
Weight of case 1.35 kg  
EAN of inner box 4038917132565  
EAN case 4038917030472

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

