



Product Specification
DishRack, height: 370 mm, green, for 88 petri dishes with 92 mm Ø

Page 1



Product description

Order number	93.1646.004
Product description	DishRack, height: 370 mm, material: PC, green, for 88 petri dishes with 92 mm Ø, 1 piece(s)/case

Product characteristics

Storage capacity	88
------------------	----

Size

Height of product	370 mm
-------------------	--------

Material & colours

Product material	Polycarbonate (PC)
Colour of product	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**DishRack, height: 370 mm, green, for 88 petri
dishes with 92 mm Ø**

Page 2

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest case
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 1
Piece(s) / pallet 108
Depth of case 383 mm
Width of case 193 mm
Height of case 190 mm
Case volume 0.014 cbm
Weight of product 1.086 kg
Weight of case 1.35 kg
EAN of inner box 4038917132589
EAN case 4038917030496

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT