



## Product description

Order number	93.1647.001
Product description	DishRack, height: 240 mm, material: PC, red, for 52 petri dishes with 92 mm Ø, 1 piece(s)/case

## Product characteristics

Storage capacity	52
------------------	----

### Size

Height of product	240 mm
-------------------	--------

## Material & colours

Product material	Polycarbonate (PC)
Colour of product	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**DishRack, height: 240 mm, red, for 52 petri dishes  
with 92 mm Ø**

**Purity & certification**

Product category                   no medical device | IVD  
Batched                               no

**Packaging**

Minimum order qty.               1  
Type of smallest                   case  
subpackaging  
Piece(s) / inner box               1  
Piece (s) / outer case           1  
Piece(s) / pallet                 162  
Depth of case                     253 mm  
Width of case                     193 mm  
Height of case                    190 mm  
Case volume                       0.0093 cbm  
Weight of product                0.617 kg  
Weight of case                    0.8 kg  
EAN of inner box                 4038917132602  
EAN case                         4038917030519

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

