



### **Product description**

Order number	93.1647.050
Product description	Silicone strip, for DishRack 50, length: 240 mm, material: silicone, translucent, 1 piece(s)/bag

### **Product characteristics**

Suitable for	DishRack 50
--------------	-------------

### **Size**

Width of product	65 mm
Length of product	240 mm

### **Material & colours**

Product material	silicone
Colour of product	translucent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

### **Packaging**

Minimum order qty.                      1  
Type of smallest                              bag  
subpackaging  
Piece(s) / inner box                      1  
Piece (s) / outer case                      30  
Piece(s) / pallet                              4860  
Depth of case                              253 mm  
Width of case                                193 mm  
Height of case                               190 mm  
Case volume                                0.0093 cbm  
Weight of product                        0.0606 kg  
Weight of case                              2 kg  
EAN of inner box                        4038917132640  
EAN case                                    4038917040648

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

