



**Product Specification**  
**Semi-micro cuvette, 3 ml, (HxW): 45 x 12 mm,**  
**PMMA, transparent, optical sides: 2**

**Page 1**



### Product description

Order number	67.740
Product description	Semi-micro cuvette, max. working volume: 3 ml, material: PMMA, transparent, (HxW): 45 x 12 mm, optical sides: 2, path length: 10 mm, sorted by mould cavity number, 100 piece(s)/polystyrene container

### Product characteristics

Volume of work	1 ml - 3 ml
Number of optical sides	2
Path length	10 mm
Wavelength in the cuvette area	300 nm

### Size

Max. volume	3 ml
Width of product	12 mm
Height of product	45 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Semi-micro cuvette, 3 ml, (HxW): 45 x 12 mm,**  
**PMMA, transparent, optical sides: 2**

**Page 2**

**Material & colours**

Product material Polymethyl methacrylate (PMMA)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Batched no

**Packaging**

Minimum order qty. 2000  
Type of smallest subpackaging polystyrene container  
Piece(s) / inner box 100  
Piece (s) / outer case 2000  
Piece(s) / pallet 48000  
Depth of case 458 mm  
Width of case 288 mm  
Height of case 415 mm  
Case volume 0.0547 cbm  
Weight of product 0.0026 kg  
Weight of case 6.4 kg  
EAN of inner box 4038917083867  
EAN case 4038917002745

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-09-01

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**