



Product Specification
Cuvette, 4.2 ml, (HxW): 45 x 12 mm, PS,
transparent, optical sides: 2

Page 1



Product description

Order number	67.745
Product description	Cuvette, max. working volume: 4.2 ml, material: PS, transparent, (HxW): 45 x 12 mm, optical sides: 2, suitable for use beginning at 330 nm, path length: 10 mm, 500 piece(s)/bag

Product characteristics

Application	suitable for use beginning at 330 nm
Volume of work	2,5 ml - 4,2 ml
Number of optical sides	2
Path length	10 mm
Wavelength in the cuvette area	330 nm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Cuvette, 4.2 ml, (HxW): 45 x 12 mm, PS,
transparent, optical sides: 2

Page 2

Size

Max. volume	4.5 ml
Width of product	12 mm
Height of product	45 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	84000
Depth of case	388 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0367 cbm
Weight of product	0.002 kg
Weight of case	4.42 kg
EAN of inner box	4038917083911
EAN case	4038917002790

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT