



Product description

| | |
|---------------------|---|
| Order number | 67.746 |
| Product description | Semi-micro cuvette, max. working volume: 3 ml, material: PS, transparent, (HxW): 45 x 12 mm, optical sides: 2, path length: 10 mm, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------------------|-------------|
| Volume of work | 1 ml - 3 ml |
| Number of optical sides | 2 |
| Path length | 10 mm |
| Wavelength in the cuvette area | 330 nm |

Size

| | |
|-------------------|-------|
| Max. volume | 3 ml |
| Width of product | 12 mm |
| Height of product | 45 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Semi-micro cuvette, 3 ml, (HxW): 45 x 12 mm, PS,
transparent, optical sides: 2

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 2000
Type of smallest
subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 2000
Piece(s) / pallet 96000
Depth of case 388 mm
Width of case 378 mm
Height of case 215 mm
Case volume 0.0315 cbm
Weight of product 0.0023 kg
Weight of case 5.02 kg
EAN of inner box 4038917083928
EAN case 4038917002806

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT