



Product description

Order number	82.1184.500
Product description	Petri dish, (ØxH): 150 x 20 mm, material: PS, transparent, with ventilation cams, sterile, 10 piece(s)/tube bag

Product characteristics

Minimum temperature	0 °C
Maximum temperature	80 °C
Ventilation cams	yes

Size

Max. volume	275 ml
Diameter	150 mm
Height of product	20 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Petri dish, 150 x 20 mm, transparent, with ventilation cams

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Latex-free yes
Purity standard sterile
Sterilisation Gamma irradiation
Safety data sheet no
Batched yes

Packaging

Minimum order qty. 100
Type of smallest subpackaging tube bag
Piece(s) / inner box 10
Piece (s) / outer case 100
Piece(s) / pallet 1800
Depth of case 438 mm
Width of case 324 mm
Height of case 540 mm
Case volume 0.0766 cbm
Weight of product 0.0524 kg
Weight of case 6.81 kg
EAN of inner box 4038917347860
EAN case 4038917347853

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT