



Product description

Order number	82.1473.040
Product description	Petri dish, (ØxH): 92 x 16 mm, material: PS, yellow, with ventilation cams, 20 piece(s)/tube bag

Product characteristics

Minimum temperature	0 °C
Maximum temperature	80 °C
Ventilation cams	yes

Size

Max. volume	80 ml
Diameter	92 mm
Height of product	16 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product yellow

Purity & certification

Product category no medical device | IVD
Latex-free yes
Safety data sheet no
Batched yes

Packaging

Minimum order qty. 480
Type of smallest tube bag
subpackaging
Piece(s) / inner box 480
Piece (s) / outer case 480
Piece(s) / pallet 9600
Depth of case 584 mm
Width of case 387 mm
Height of case 360 mm
Case volume 0.0814 cbm
Weight of product 0.0169 kg
Weight of case 9.25 kg
EAN of inner box 4038917129374
EAN case 4038917022156

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

