



### **Product description**

Order number	92.180.616
Product description	Block rotor, for 14 tubes up to 28 mm diameter, for SARMIX® M 2000, 1 piece(s)/case

### **Product characteristics**

Suitable for	SARMIX® M 2000
Storage capacity	14

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	10
Piece(s) / pallet	10

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Block rotor, for 14 tubes up to 28 mm diameter,**  
**for SARMIX® M 2000**

**Page 2**

Case volume	0.0122 cbm
Weight of product	0.18 kg
Weight of case	3.01 kg
EAN of inner box	4038917160858
EAN case	4038917160841

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**