



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 92.181.610   |
| Product description | Rotor, for 36 collection tubes, for micro tube Ø 10.8 mm |

### **Product characteristics**

|                  |                      |
|------------------|----------------------|
| Suitable for     | micro tube Ø 10.8 mm |
| Storage capacity | 36                   |

### **Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

### **Packaging**

|                        |            |
|------------------------|------------|
| Minimum order qty.     | 1          |
| Piece(s) / inner box   | 1          |
| Piece (s) / outer case | 1          |
| Piece(s) / pallet      | 40         |
| Case volume            | 0.0056 cbm |
| Weight of product      | 1.3595 kg  |
| Weight of case         | 1.36 kg    |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Rotor, for 36 collection tubes, for micro tube Ø**  
**10.8 mm**

**Page 2**

EAN of inner box                      4038917156233  
EAN case                                    4038917156226

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**