



## Product Specification

**IsoFreeze® PCR Rack, PP, format: 12 x 8, suitable for 0.1 ml and 0.2 ml PCR tubes, PCR strips and PCR plates**

**Page 1**



### Product description

|                     |                                                                                                                                                                                                    |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 95.984                                                                                                                                                                                             |
| Product description | IsoFreeze® PCR Rack, material: PP, multi-coloured, format: 12 x 8, suitable for 0.1 ml and 0.2 ml PCR tubes, PCR strips and PCR plates, Colour change from purple to pink at 7 °C, 2 piece(s)/case |

### Product characteristics

|              |                                                        |
|--------------|--------------------------------------------------------|
| Suitable for | 0.1 ml and 0.2 ml PCR tubes, PCR strips and PCR plates |
| Format       | 12 x 8                                                 |

### Material & colours

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | multi-coloured     |
| Colour of cap     | transparent        |
| Colour of frame   | white              |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

IsoFreeze® PCR Rack, PP, format: 12 x 8, suitable for 0.1 ml and 0.2 ml PCR tubes, PCR strips and PCR plates

Purity & certification

Product category                      no medical device | IVD  
Batched                                      no

Packaging

Minimum order qty.                      2  
Type of smallest subpackaging                      case  
Piece(s) / inner box                      2  
Piece (s) / outer case                      10  
Piece(s) / pallet                      1280  
Depth of case                      197 mm  
Width of case                      175 mm  
Height of case                      584 mm  
Case volume                      0.0201 cbm  
Weight of product                      0.2144 kg  
Weight of case                      2.95 kg  
EAN of inner box                      4038917373425  
EAN case                      4038917373418

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

