



### **Product description**

Order number	95.988.002
Product description	Base station, material: PC, blue, (LxWxH): 128 x 85.5 x 30.4 mm, suitable for work tray 95.987.002 & 0.2 ml PCR tubes, 5 piece(s)/case

### **Product characteristics**

Suitable for	work tray 95.987.002 & 0.2 ml PCR tubes
--------------	---

### **Size**

Width of product	85.5 mm
Height of product	30.4 mm
Length of product	128 mm

### **Material & colours**

Product material	Polycarbonate (PC)
Colour of product	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      5  
Type of smallest  
subpackaging                              case  
Piece(s) / inner box                      5  
Piece (s) / outer case                      75  
Piece(s) / pallet                            3600  
Depth of case                              598 mm  
Width of case                                398 mm  
Height of case                               145 mm  
Case volume                                0.0345 cbm  
Weight of product                        0.0302 kg  
Weight of case                              2.83 kg  
EAN of inner box                         4038917124171  
EAN case                                      4038917019132

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

