



Product description

Order number	72.991.102
Product description	PCR strip of 8, High Profile, PCR Performance Tested, transparent, material: PP, flat cap, cap attached, 120 piece(s)/bag

Product characteristics

Graduation	no
Cap	flat cap
Closure type	snap-on lid
Base shape	conical
High or low profile	High Profile
Number of wells	8
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	20
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
PCR strip of 8, PCR Performance Tested,
transparent, PP, flat cap

Page 2

Size

Sample volume 0.2 ml

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard PCR Performance Tested
Batched yes

Packaging

Minimum order qty. 1200
Type of smallest bag
subpackaging
Piece(s) / inner box 120
Piece (s) / outer case 1200
Piece(s) / pallet 40320
Case volume 0.0354 cbm
Weight of product 0.0018 kg
Weight of case 2.77 kg
EAN of inner box 4038917231657
EAN case 4038917231640

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT