



### Product description

Order number	72.991.103
Product description	PCR strip of 8, up to 100 µl, Low Profile, PCR Performance Tested, transparent, material: PP, flat cap, cap attached, suitable for qPCR, 12 piece(s)/minigrip bag

### Product characteristics

Graduation	no
Cap	flat cap
Closure type	push cap
Suitable for qPCR	yes
Base shape	conical
High or low profile	Low Profile
Number of wells	8
Volume of work	5 µl - 0,1 ml
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Sample volume	0.1 ml
Max. volume	100 µl
Diameter	6.26 mm
Width of product	10.6 mm
Height of product	15.4 mm
Length of product	71.6 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	PCR Performance Tested
Batched	yes

### **Packaging**

Minimum order qty.	120
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	120
Piece (s) / outer case	1200
Piece(s) / pallet	43200
Depth of case	305 mm
Width of case	279 mm
Height of case	305 mm
Case volume	0.026 cbm
Weight of product	0.001 kg
Weight of case	2.2 kg
EAN of inner box	4038917231688
EAN case	4038917231671

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

