



Product description

Order number	95.1993
Product description	PCR film, free from DNase/RNase, material: PL, transparent, suitable for qPCR, self-adhesive, application area: -40 - 120 °C, thin material, high optical clarity, 100 piece(s)/bag

Product characteristics

Closure type	film
Light protection	no
Suitable for qPCR	yes
Minimum temperature	-40 °C
Maximum temperature	120 °C

Size

Width of product	80 mm
Length of product	137 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material polyester (PL)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard free from DNase/RNase
Batched yes

Packaging

Minimum order qty. 100
Type of smallest
subpackaging bag
Piece(s) / inner box 100
Piece (s) / outer case 1600
Piece(s) / pallet 192000
Depth of case 216 mm
Width of case 216 mm
Height of case 203 mm
Case volume 0.0095 cbm
Weight of product 0.002 kg
Weight of case 3.74 kg
EAN of inner box 4038917303569
EAN case 4038917388221

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

