



Product description

Order number	93.838
Product description	Rack, material: PC, transparent, format: 12 x 4, (LxWxH): 257 x 90 x 40 mm, for 48 collection tubes, suitable for tubes, S-Monovette® 11 mm Ø, 1 piece(s)/case

Product characteristics

Suitable for	tubes, S-Monovette® 11 mm Ø
Storage capacity	48
Format	12 x 4
Autoclavable (max. °C)	121 °C

Size

Width of product	90 mm
Height of product	40 mm
Length of product	257 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Rack, PC, format: 12 x 4, suitable for tubes,
S-Monovette® 11 mm Ø

Page 2

Material & colours

Product material Polycarbonate (PC)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest case
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 40
Piece(s) / pallet 1200
Depth of box 92 mm
Width of box 42 mm
Height of box 259 mm
Depth of case 398 mm
Width of case 398 mm
Height of case 340 mm
Case volume 0.0539 cbm
Weight of product 0.2138 kg
Weight of case 10.29 kg
EAN of inner box 4038917084352
EAN case 4038917003216

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT