



Product description

Order number	93.844.100
Product description	Rack, material: PC, transparent, format: 10 x 2, (LxWxH): 257 x 62 x 55 mm, for 20 collection tubes, suitable for tubes, all S-Monovette® diameters, 1 piece(s)/case

Product characteristics

Suitable for	tubes, all S-Monovette® diameters
Storage capacity	20
Format	10 x 2
Autoclavable (max. °C)	121 °C

Size

Width of product	62 mm
Height of product	55 mm
Length of product	257 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polycarbonate (PC)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest case
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 42
Piece(s) / pallet 1260
Depth of box 65 mm
Width of box 58 mm
Height of box 260 mm
Depth of case 398 mm
Width of case 398 mm
Height of case 340 mm
Case volume 0.0539 cbm
Weight of product 0.1907 kg
Weight of case 10.14 kg
EAN of inner box 4038917121712
EAN case 4038917017961

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

