



Product description

Order number	95.1439
Product description	Rack, material: styrofoam, white, format: 20 x 5, (LxWxH): 338 x 90 x 50 mm, suitable for S-Monovette® Ø 11 mm, 1 piece(s)/blister

Product characteristics

Suitable for	S-Monovette® Ø 11 mm
Format	20 x 5

Size

Width of product	90 mm
Height of product	50 mm
Length of product	338 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Rack, styrofoam, format: 20 x 5, suitable for
S-Monovette® Ø 11 mm

Page 2

Material & colours

Product material styrofoam
Colour of product white

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest
subpackaging blister
Piece(s) / inner box 1
Piece (s) / outer case 40
Piece(s) / pallet 640
Depth of case 498 mm
Width of case 428 mm
Height of case 380 mm
Case volume 0.081 cbm
Weight of product 0.0235 kg
Weight of case 1.7 kg
EAN of inner box 4038917086646
EAN case 4038917005173

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT