



### **Product description**

Order number	95.64.213
Product description	Rack, material: styrofoam, white, format: 20 x 5, (LxWxH): 336 x 89 x 49 mm, for 100 collection tubes, suitable for tubes Ø 13 mm, 1 piece(s)/blister

### **Product characteristics**

Suitable for	tubes Ø 13 mm
Storage capacity	100
Format	20 x 5

### **Size**

Width of product	89 mm
Height of product	49 mm
Length of product	336 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	styrofoam
Colour of product	white
Closure material	styrofoam
Colour of cap	white

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	35
Piece(s) / pallet	700
Depth of case	458 mm
Width of case	408 mm
Height of case	350 mm
Case volume	0.0654 cbm
Weight of product	0.03 kg
Weight of case	1.75 kg
EAN of inner box	4038917095563
EAN case	4038917030809

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

