



### **Product description**

|                     |                                                                                                                                                               |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 95.64.250                                                                                                                                                     |
| Product description | Rack, material: styrofoam, white, format: 20 x 5, (LxWxH): 416 x 107 x 50 mm, for 100 collection tubes, suitable for tubes Ø 15 and 16 mm, 1 piece(s)/blister |

### **Product characteristics**

|                  |                      |
|------------------|----------------------|
| Suitable for     | tubes Ø 15 and 16 mm |
| Storage capacity | 100                  |
| Format           | 20 x 5               |

### **Size**

|                   |        |
|-------------------|--------|
| Width of product  | 107 mm |
| Height of product | 50 mm  |
| Length of product | 416 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

|                   |           |
|-------------------|-----------|
| Product material  | styrofoam |
| Colour of product | white     |
| Closure material  | styrofoam |
| Colour of cap     | white     |

### **Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | no                      |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1             |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 1             |
| Piece (s) / outer case        | 28            |
| Piece(s) / pallet             | 448           |
| Depth of case                 | 498 mm        |
| Width of case                 | 428 mm        |
| Height of case                | 380 mm        |
| Case volume                   | 0.081 cbm     |
| Weight of product             | 0.0309 kg     |
| Weight of case                | 1.65 kg       |
| EAN of inner box              | 4038917095587 |
| EAN case                      | 4038917007764 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

