



### **Product description**

Order number	95.983
Product description	IsoFreeze® MCT Rack, material: PP, multi-coloured, format: 6 x 4, suitable for reaction tubes Ø 10.8 mm (1.5 ml & 2 ml), 1 piece(s)/case

### **Product characteristics**

Label/ Print	with print
Suitable for	reaction tubes Ø 10.8 mm (1.5 ml & 2 ml)
Format	6 x 4

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	multi-coloured
Colour of cap	transparent
Colour of frame	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category no medical device | IVD  
Batched no

### **Packaging**

Minimum order qty. 1  
Type of smallest case  
subpackaging  
Piece(s) / inner box 1  
Piece (s) / outer case 5  
Piece(s) / pallet 240  
Depth of case 197 mm  
Width of case 93 mm  
Height of case 584 mm  
Case volume 0.0107 cbm  
Weight of product 0.349 kg  
Weight of case 2.5 kg  
EAN of inner box 4038917394895  
EAN case 4038917394888

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

