



Product description

| | |
|---------------------|---|
| Order number | 95.64.924 |
| Product description | Storage box, slip-on lid, material: cardboard, white, format: 12 x 12, for 144 collection tubes, suitable for tubes up to max. 9.5 mm Ø, tube height between 19-26 mm, 1 piece(s)/bag |

Product characteristics

| | |
|------------------|---|
| Closure type | slip-on lid |
| Suitable for | tubes up to max. 9.5 mm Ø, tube height between 19-26 mm |
| Storage capacity | 144 |
| Format | 12 x 12 |

Size

| | |
|-------------------|--------|
| Width of product | 135 mm |
| Height of product | 28 mm |
| Length of product | 135 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



Product Specification
Storage box, slip-on lid, cardboard, format: 12 x 12, for 144 collection tubes

Page 2

Material & colours

| | |
|-------------------|-----------|
| Product material | cardboard |
| Colour of product | white |
| Closure material | cardboard |
| Colour of cap | white |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 216 |
| Piece(s) / pallet | 2592 |
| Depth of case | 521 mm |
| Width of case | 521 mm |
| Height of case | 610 mm |
| Case volume | 0.1653 cbm |
| Weight of product | 0.0632 kg |
| Weight of case | 15.6 kg |
| EAN of inner box | 4038917095723 |
| EAN case | 4038917067461 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT