



Product description

Order number 95.64.998
Product description Storage box, slip-on lid, material: cardboard, white

Product characteristics

Label/ Print with print
Closure type slip-on lid

Size

Width of product 135 mm
Height of product 45 mm
Length of product 135 mm

Material & colours

Product material cardboard
Colour of product white
Closure material cardboard
Colour of cap white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 10
Piece(s) / inner box 10
Piece (s) / outer case 120
Piece(s) / pallet 2160
Case volume 0.1134 cbm
Weight of product 0.065 kg
Weight of case 7.86 kg
EAN of inner box 4038917246927
EAN case 4038917246910

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

