



### **Product description**

Order number	65.1379.007
Product description	Cap, violet, suitable for tubes Ø 10-17 mm, 1,000 piece(s)/bag

### **Product characteristics**

Closure type	push cap
Suitable for	tubes Ø 10-17 mm

### **Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	violet

### **Purity & certification**

Product category	no medical device   IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	360000
Case volume	0.0492 cbm
Weight of product	0 kg
Weight of case	5.35 kg
EAN of inner box	4038917128025
EAN case	4038917021555

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

