



### **Product description**

Order number	65.803
Product description	Push cap, natural, suitable for tubes Ø 15.7 mm & cuvette 67.743, 1,000 piece(s)/bag

### **Product characteristics**

Cap	standard cap
Closure type	push cap
Suitable for	tubes Ø 15.7 mm & cuvette 67.743

### **Size**

Diameter	16.35 mm
Height of product	20.45 mm

### **Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category no medical device | IVD  
Batched no

**Packaging**

Minimum order qty. 5000  
Type of smallest subpackaging bag  
Piece(s) / inner box 1000  
Piece (s) / outer case 5000  
Piece(s) / pallet 210000  
Depth of case 398 mm  
Width of case 398 mm  
Height of case 250 mm  
Case volume 0.0396 cbm  
Weight of product 0.0013 kg  
Weight of case 7.1 kg  
EAN of inner box 4038917084192  
EAN case 4038917003056

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

