



### **Product description**

Order number	65.806
Product description	Push cap, natural, suitable for tubes Ø 13 mm, 1,000 piece(s)/bag

### **Product characteristics**

Cap	standard cap
Closure type	push cap
Suitable for	tubes Ø 13 mm

### **Size**

Diameter	13.3 mm
Height of product	15.2 mm

### **Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

### **Packaging**

Minimum order qty.                      10000  
Type of smallest  
subpackaging                              bag  
Piece(s) / inner box                      1000  
Piece (s) / outer case                      10000  
Piece(s) / pallet                            420000  
Depth of case                              398 mm  
Width of case                                398 mm  
Height of case                               250 mm  
Case volume                                0.0396 cbm  
Weight of product                        0.0007 kg  
Weight of case                              7.36 kg  
EAN of inner box                          4038917084222  
EAN case                                    4038917003087

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

