



Product description

| | |
|---------------------|---|
| Order number | 55.1579.002 |
| Product description | Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, suitable for flow cytometer like FACSCanto, FACS Aria, FACSCalibur, LRS2, FC500, sterile, 50 piece(s)/bag |

Product characteristics

| | |
|--------------------------|--|
| Graduation | no |
| Closure type | push cap |
| Suitable for | flow cytometer like FACSCanto, FACS Aria, FACSCalibur, LRS2, FC500 |
| Base shape | round base |
| Volume of work | 5 ml |
| Centrifugation max (RCF) | 4000 x g |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|-------|
| Diameter | 12 mm |
| Length of product | 75 mm |

Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | sealed at 0.6 bar/10 min |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 56000 |
| Depth of case | 598 mm |
| Width of case | 398 mm |
| Height of case | 125 mm |
| Case volume | 0.0298 cbm |
| Weight of product | 0.0026 kg |
| Weight of case | 3.18 kg |
| EAN of inner box | 4038917162746 |
| EAN case | 4038917162739 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

