



### **Product description**

Order number	55.467.300
Product description	Tube, working volume: 8 ml, (LxØ): 100 x 13 mm, material: PS, round base, transparent, push cap, 500 piece(s)/StackPack

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	8 ml
Centrifugation max (RCF)	3500 x g

### **Size**

Diameter	13 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material                      Polystyrene (PS)  
Colour of product                      transparent

### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      yes

### **Packaging**

Minimum order qty.                      500  
Type of smallest  
subpackaging                              StackPack  
Piece(s) / inner box                      500  
Piece (s) / outer case                      2000  
Piece(s) / pallet                              96000  
Case volume                                  0.0366 cbm  
Weight of product                          0.0039 kg  
Weight of case                                8.86 kg  
EAN of inner box                              4038917367042  
EAN case                                        4038917367035

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

