



### **Product description**

Order number	55.468.001
Product description	Tube, working volume: 13 ml, (LxØ): 95 x 16.8 mm, material: PS, round base, transparent, push cap, natural, cap assembled, sterile, 100 piece(s)/bag

### **Product characteristics**

Graduation	no
Cap	standard cap
Closure type	push cap
Base shape	round base
Volume of work	13 ml
Centrifugation max (RCF)	4000 x g

### **Size**

Diameter	16.8 mm
Length of product	95 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Colour of cap	natural

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### **Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	28000
Depth of case	598 mm
Width of case	398 mm
Height of case	125 mm
Case volume	0.0298 cbm
Weight of product	0.0069 kg
Weight of case	4.03 kg
EAN of inner box	4038917106962
EAN case	4038917011594

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

