



Product description

Order number	55.516.015
Product description	Tube, working volume: 8 ml, (LxØ): 100 x 13 mm, material: PP, round base, transparent, push cap, 500 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	8 ml
Centrifugation max (RCF)	3500 x g
Autoclavable (max. °C)	121 °C

Size

Diameter	13 mm
Height of product	100 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	StackPack
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	96000
Case volume	0.0366 cbm
Weight of product	0.0041 kg
Weight of case	9.25 kg
EAN of inner box	4038917108959
EAN case	4038917074681

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

